



- Ergonomic design, comfortable pipetting
- Wide volume range (0.1-10000 μ l)
- Single pipettes above 10 μ l are equipped with replaceable hydrophobic cores
- Tips ejection stick is easy to be disassembled and autoclavable
- Lightweight, durable and precise

Technology Creates The Future

Genex pipette

BIO-DL Genex adjustable pipette

Order Information



Item No.	Channels	Range	Increment	Measuring Volume	Inaccuracy	Imprecision	Compatible Tip
72441011	1	0.1-2.5µl	0.05µl	2.5µl	±2.50%	2.00%	10µl
				1.25µl	±3.00%	3.00%	
				0.25µl	±12.00%	6.00%	
72441031	1	0.5-10µl	0.1µl	10µl	±1.00%	0.80%	10µl
				5µl	±1.50%	1.50%	
				1µl	±2.50%	1.50%	
72441041	1	2-20µl	0.5µl	20µl	±0.90%	0.40%	200µl & 300µl
				10µl	±1.20%	1.00%	
				2µl	±3.00%	2.00%	
72441051	1	5-50µl	0.5µl	50µl	±0.60%	0.30%	200µl & 300µl
				25µl	±0.90%	0.60%	
				5µl	±2.00%	2.00%	
72441061	1	10-100µl	1µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72441071	1	20-200µl	1µl	200µl	±0.60%	0.20%	200µl & 300µl
				100µl	±0.80%	0.30%	
				20µl	±2.50%	0.80%	
72441081	1	100-1000µl	5µl	1000µl	±0.60%	0.20%	1000µl
				500µl	±0.70%	0.25%	
				100µl	±2.00%	0.70%	
72441101	1	1000-5000µl	50µl	5000µl	±0.50%	0.20%	5000µl
				2500µl	±0.60%	0.30%	
				1000µl	±0.70%	0.30%	
72441111	1	1-10ml	100µl	10000µl	±0.60%	0.20%	10ml
				5000µl	±1.20%	0.30%	
				1000µl	±3.00%	0.60%	



Item No.	Channels	Range	Increment	Measuring Volume	Inaccuracy	Imprecision	Compatible Tip
72411031	8	0.5-10µl	0.1µl	10µl	±1.50%	1.50%	10µl
				5µl	±2.50%	2.50%	
				1µl	±4.00%	4.00%	
72411051	8	5-50µl	0.5µl	50µl	±1.00%	0.50%	200µl & 300µl
				25µl	±1.50%	1.00%	
				5µl	±3.00%	2.00%	
72411061	8	10-100µl	0.5µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72411141	8	50-300µl	5µl	300µl	±0.70%	0.25%	300µl & 350µl
				150µl	±1.00%	0.50%	
				50µl	±1.50%	0.80%	
72411081	8	100-1000µl	10µl	1000µl	±0.60%	0.20%	1000µl
				500µl	±0.70%	0.25%	
				100µl	±2.00%	0.70%	
72421031	12	0.5-10µl	0.1µl	10µl	±1.50%	1.50%	10µl
				5µl	±2.50%	2.50%	
				1µl	±4.00%	4.00%	
72421051	12	5-50µl	0.5µl	50µl	±1.00%	0.50%	200µl & 300µl
				25µl	±1.50%	1.00%	
				5µl	±3.00%	2.00%	
72421061	12	10-100µl	0.5µl	100µl	±0.80%	0.20%	200µl & 300µl
				50µl	±1.00%	0.40%	
				10µl	±3.00%	1.00%	
72421141	12	50-300µl	5µl	300µl	±0.70%	0.25%	300µl & 350µl
				150µl	±1.00%	0.50%	
				50µl	±1.50%	0.80%	

BIO-DL Genex fixed single channel pipette

Item No.	Channels	Range	Inaccuracy	Imprecision	Compatible Tip
72431211	1-ch	5µl	±0.30%	1.2%	10µl
72431221	1-ch	10µl	±0.80%	0.8%	10µl
72431231	1-ch	20µl	±0.50%	0.5%	200µl
72431241	1-ch	25µl	±0.30%	0.3%	200µl
72431251	1-ch	50µl	±0.30%	0.3%	200µl
72431261	1-ch	100µl	±0.30%	0.3%	200µl
72431271	1-ch	200µl	±0.20%	0.2%	200µl
72431281	1-ch	250µl	±0.20%	0.2%	1000µl
72431291	1-ch	500µl	±0.20%	0.2%	1000µl
72431311	1-ch	1000µl	±0.20%	0.2%	1000µl
72431321	1-ch	2000µl	±0.15%	0.15%	5000µl
72431331	1-ch	5000µl	±0.15%	0.15%	5000µl

